**REF**

160245, 151972, 161175, 161188

English | 624591/07  
2019-09-30

---

	Copyright and conventions .....	2
1	Safety Information .....	2
2	Overview .....	2
3	Reprocessing the expiratory valve set .....	3
4	Using a reprocessed expiratory valve set .....	5
5	Disposal .....	5
6	Parts and accessories .....	5

## Copyright and conventions

© 2019 Hamilton Medical AG. All rights reserved. Printed in Switzerland.

No part of this publication may be reproduced, stored in a database or retrieval system, or transmitted in any form or by any means, electronic, mechanical, or by photocopying, recording, or otherwise, without prior written permission of Hamilton Medical AG.

This document may be revised, replaced, or made obsolete by other documents by Hamilton Medical AG at any time and without notice. Ensure that you have the most current applicable version of this document; if in doubt, contact the technical support department of Hamilton Medical AG, Switzerland. While the information set forth herein is believed to be accurate, it is not a substitute for the exercise of professional judgment.

Nothing in this document shall limit or restrict in any way Hamilton Medical AG's right to revise or otherwise change or modify the equipment (including its software) described herein, without notice. In the absence of an express, written agreement to the contrary, Hamilton Medical AG has no obligation to furnish any such revisions, changes, or modifications to the owner or user of the equipment (including software) described herein.

The equipment must be operated and serviced only by trained professionals.

For all proprietary as well as third-party trademarks used by Hamilton Medical AG, see [www.hamilton-medical.com/trademarks](http://www.hamilton-medical.com/trademarks). Product and/or company names marked with a <sup>§</sup> symbol may be the trademarks and/or registered trademarks of their respective owners, including but not limited to Aerogen<sup>§</sup>, Nihon Kohden<sup>§</sup>, Masimo<sup>§</sup>, Masimo SET<sup>§</sup>, Masimo rainbow SET<sup>§</sup>, and Capnostat<sup>§</sup>.

### Conventions used in this guide

#### **WARNING**

Alerts the user to the possibility of injury, death, or other serious adverse reactions associated with the use or misuse of the device.

#### **CAUTION**

Alerts the user to the possibility of a problem with the device associated with its use or misuse, such as device malfunction, device failure, damage to the device, or damage to other property.

#### **NOTICE**

Emphasizes information of particular importance.

## 1 Safety Information

#### **WARNING**

- Clean and sterilize the component directly after use.
- Hamilton Medical cannot be held responsible for the correct functioning of components that are not reprocessed according to these instructions.
- Ensure that only processes that have been specifically validated for the product or device are used, and that the validated parameters are used with every cycle.
- Any used component must be handled as a contaminated item. Follow all local, state, and federal regulations with respect to environmental protection when disposing of used components.
- Follow hospital infection control procedures, as well as local laws and regulations. This applies in particular to the various regulations regarding an effective deactivation of prions.
- Defective components or components that fail the tightness test must not be used.
- Visually inspect the expiratory valve set components for damage prior to use.

#### **CAUTION**

- *Autoclavable components have a limited life span.*
- *If reprocessing components manually, take care that no damage is caused by the use of hard brushes, scouring agents, or by the exertion of too much force.*
- *The use of rinse aids will reduce the life span of the components, as it can lead to early failure and cracks in the plastic.*
- *(USA only) Only use EPA registered cleaning solutions.*

#### **NOTICE**

Make sure that the reprocessing does not damage the steel ring and the membrane.

## 2 Overview

The expiratory valve set must be cleaned and sterilized after each use in accordance with EN ISO 17664.


Rinse aids will cause premature damage and reduce product life span, and should not be used. Hamilton Medical does not guarantee the expiratory valve set's life span if rinse aids are used.

### 2.1 Components and materials

The following section provides an overview of the different expiratory valve sets, their components, and the materials they comprise of. This recommendation is valid for the Hamilton Medical products listed below.

All materials are heat resistant up to 140°C (284°F).


#### **HAMILTON-C1/T1/MR1 expiratory valve set (Adult/Ped shown)**

- 1 Membrane
- 2 Housing
- 3 Locking ring

PN	Component	Materials
161175 (Adult/Ped)	Housing	Polycarbonate
	Locking ring	Polyamide 12
	Membrane	Silicone rubber
161188 (Neonatal)	Metal plate on membrane	Stainless steel
	Housing	Polycarbonate
	Locking ring	Polyamide 12
	Membrane	Silicone rubber
	Metal plate on membrane	Stainless steel


### HAMILTON-C2/C3/C6 expiratory valve set

- 1 Membrane
- 2 Housing
- 3 Locking ring<sup>1</sup>

PN	Component	Materials
160245	Housing	Polyamide 12
	Locking ring	Polyamide 12
	Membrane	Silicone rubber
	Metal plate on membrane	Stainless steel

### HAMILTON-G5/S1 expiratory valve set

- 1 Membrane
- 2 Housing

PN	Component	Materials
151972	Housing	Polysulfone
	Membrane	Silicone rubber
	Metal plate on membrane	Stainless steel

## 3 Reprocessing the expiratory valve set

Reprocess the expiratory valve set immediately after use. Reprocessing the expiratory valve set comprises the following steps:

To reprocess the expiratory valve set ...	See ...
Disassemble	Section 3.1
Clean	Section 3.2
Carry out a visual test	Section 3.3
Package	Section 3.4
Sterilize	Section 3.5
Store for future use	Section 3.6


### 3.1 Disassembly

#### To disassemble the components


- (HAMILTON-C1/T1/MR1 only) Remove the locking ring from the housing by twisting the ring counter-clockwise (1) and pulling it off of the housing (2).
- Holding the expiratory valve housing, remove the silicone membrane by lifting it up (3).

Figure 1. Disassembling the expiratory valve set


HAMILTON-C1/T1/MR1 (Adult/Ped shown)

HAMILTON-C2/C3/C6

HAMILTON-G5/S1

### 3.2 Cleaning

Remove macroscopic impurities from the expiratory valve set by rinsing or wiping. You must not use any hard tools or hard brushes to remove resilient impurities.

You can clean the expiratory valve using either of the following methods:

To perform ...	See ...
Manual cleaning	Section 3.2.1
Mechanical cleaning	Section 3.2.2

#### 3.2.1 Manual cleaning

The following table lists the recommended cleaning agents to use for manual cleaning of the expiratory valve set.

Table 1. Recommended cleaning agent, manual cleaning

Manufacturer	Product	Concentration
Dr. Weigert	Neodisher Mediclean forte <sup>5</sup>	2%

For details about the cleaning agent, refer to the manufacturer's *Instructions for use*.

#### To manually clean the components

- Disassemble the expiratory valve set (Section 3.1).
- Rinse the expiratory valve set components at least five times with the recommended cleaning agent (Table 1).
- Soak the expiratory valve set components in the recommended cleaning solution (Table 1) for 30 minutes at a temperature of 40°C (104°F).  
Make sure that all parts of the expiratory valve set components are fully submerged in the solution.
- Rinse the expiratory valve set components at least five times with the recommended cleaning agent (Table 1).
- Remove matter and larger exterior impurities by carefully scrubbing the expiratory valve set components with a soft brush or soft towel.
- Rinse the expiratory valve components at least five times in deionized water.
- Repeat the cleaning process if the last cleaning solution was not clear or there are still visible impurities on the expiratory valve components.

<sup>1</sup> The locking ring for the HAMILTON-C2/C3/C6 expiratory valve set is permanently attached to the housing.

### 3.2.2 Mechanical cleaning

The expiratory valve set must be reprocessed in such a manner that hygienic and safe reuse can be assured. Mechanical cleaning should only be carried out in a device that complies with ISO 15883<sup>2</sup>. Place the expiratory valve set components in such a manner that they are easy to clean and the effectiveness of cleaning is not impaired.

Hamilton Medical recommends the DES-VAR-TD-Anaesthesia program. The following table lists the recommended device and cleaning agents to use for mechanical cleaning of the expiratory valve set.

Manufacturer	Product	Concentration
<b>Device</b>		
Miele	PG8536 Washer-disinfector	n/a
	E436/6 Injector tray	n/a
<b>Cleaning agents</b>		
Dr. Weigert	Neodisher Mediclean forte <sup>5</sup>	1%
<b>Neutralizer</b>		
Dr. Weigert	Neodisher Z <sup>5</sup>	0.1%


For details about the equipment or cleaning agents, see the manufacturer's *Instructions for use*.

Figure 2. Arranging the expiratory valve set components in the device


HAMILTON-C1/T1/MR1 (Adult/Ped shown)

HAMILTON-C2/C3/C6

HAMILTON-G5/S1

Follow the instructions provided in the *Operator's Manual* of the cleaning and disinfection device. To ensure safe cleaning, the expiratory valve housing must be connected to the corresponding receptors. The expiratory valve housing must not disconnect from the receptor during reprocessing.

Expiratory valve set components that disconnect during reprocessing must be processed again. After the cleaning process is complete, check that the expiratory valve set is completely dry and undamaged. Damaged expiratory valve sets must be discarded.

The following program parameters must be met for successful mechanical cleaning:

Table 2. Mechanical cleaning program parameters

Program parameters	Specification
Pre-rinse	One cycle using cold water for 1 min
Cleaning	One cycle at 55°C (131°F) for 5 min
Optional neutralization	One cycle using cold water for 1 min
Rinsing	One cycle using cold water for 1 min
Thermic disinfection	One cycle at 83°C (181.4°F) for 10 min
Drying	100°C (212°F) for 10 min and 95°C (203°F) for 30 min

### 3.3 Visual test

After each cleaning cycle, the expiratory valve set must be macroscopically clean, that is, free of visible residual matter and other impurities. If it is not, the entire cleaning process must be repeated.

Visually check for external damage, such as cracks, broken or deformed parts, or discoloration.

### 3.4 Packaging

Make sure that the expiratory valves are not moist during packaging.

The packaging must conform to ISO 11607 and be suitable for vapor sterilization (heat resistance up to 141°C (285.8°F)) and be sufficiently permeable to vapor.

Hamilton Medical recommends to package the components twice.

(USA only) For packaging, only use an FDA approved sterilization pouch (e.g. Medipack see through pouches).

For details on working with the packaging, see the manufacturer's *Instructions for use*.

<sup>2</sup> Not applicable in the USA.

### 3.5 Sterilization

Prior to sterilization, the expiratory valve set must be cleaned (Section 3.2).

Sterilize the expiratory valve set after cleaning and prior to reuse according to the following:

- 134.0°C (273.2°F) with prevacuum, with an exposure time of 3 minutes and a drying time of 20 minutes

Place the expiratory valve set components horizontally into the sterilizer; do not stack them.

### 3.6 Storage conditions

Storage is subject to the same guidelines as the Hamilton Medical ventilator, as specified in your ventilator *Operator's Manual*.

The expiratory valve set must be reprocessed at least every 2 years.

## 4 Using a reprocessed expiratory valve set

After each reprocessing cycle, the expiratory valve set components must be inspected for damage. If any changes are visible, the expiratory valve set must be discarded. The expiratory valve set must pass the tightness test before use (Section 4.2).


### 4.1 Reassembly

#### To reassemble the components

1. Holding the expiratory valve housing, seat the silicone membrane onto the housing (1). See Figure 3.  
The metal plate must face up and be visible.
2. (*HAMILTON-C1/T1/MR1 only*) Align the locking ring with the bottom of the housing (2) and twist clockwise until it locks into place (3).  
Be sure not to over-tighten the expiratory valve locking ring.

Figure 3. Reassembling the expiratory valve set


HAMILTON-C1/T1/MR1 (Adult/Ped shown)

HAMILTON-C2/C3/C6

HAMILTON-G5/S1

### 4.2 Testing before use

Perform a tightness test, as described in the ventilator *Operator's Manual*. If the test fails, it may be repeated once. The expiratory valve set must be replaced if the tightness test fails the second time.

### 4.3 Component life span

As long as the expiratory valve set passes the tightness test during the preoperational check, the expiratory valve set can be used. Tests and calibrations must be carried out as specified in your ventilator *Operator's Manual*.

## 5 Disposal

A used expiratory valve set must be handled as a contaminated item. Follow all local, state, and federal regulations with respect to environmental protection when disposing of used expiratory valve sets.

## 6 Parts and accessories

The following table lists the parts available for the expiratory valve set. For ordering information, refer to the e-catalog on the Hamilton Medical website or contact your Hamilton Medical representative.

Table 3. Expiratory valve set parts and accessories

PN	Description
<b>HAMILTON-C1/T1/MR1 expiratory valve set (Adult/Ped and Neonatal)</b>	
161390	Expiratory valve membrane for HAMILTON-C1/T1/MR1, box of 5
<b>HAMILTON-C2/C3/C6 expiratory valve set</b>	
160500	Expiratory valve membrane for HAMILTON-C2/C3/C6, box of 5
<b>HAMILTON-G5/S1 expiratory valve set</b>	
151233	Expiratory valve membrane for HAMILTON-G5/S1, box of 5
151228	Expiratory valve housing for HAMILTON-G5/S1, box of 1

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